

ELSODENT

Clinical and more biocompatible products

G-ETCH

High quality etching gel with an excellent contrast



Features

- Excellent balance between thixotropy and water solubility: remains in place but wash off quickly and easily.
- Non-drying, medium viscosity gel.
- Excellent contrast.



Etching of the lingual cavities.

References & presentations

JKB-50 (Blue)	"Jumbo kit" : 1 bulk pack 50 ml storage syringe of 38% phosphoric acid blue etching gel + 5 empty 3 ml delivery syringes + 40 disposable needle application tips.
JKV-50 (Green)	
JRB-50 (Blue)	Jumbo kit Refill: 1 bulk pack 50ml storage syringe of 38% phosphoric acid blue etching gel.
JRV-50 (Green)	
MET-6	"Minikit": 4 x 1,5 g syringes of 38% phosphoric acid etching gel + 10 needle tips. Only available in blue.
MET-24	24 x 1,5 g of 38% phosphoric acid etching gel. Only available in blue.
LTL	100 black needle tips (for Jumbo kits). G22.
ETL	100 blue needle tips (for Minikit). G25.

MD: class IIA. Certifying body: SNCH. Manufacturer: G-Pharma



Phone: +33 (0)1 30 37 75 75 - elsodent@elsodent.com www.elsodent.com
ELSODENT, 35 Avenue du Gros Chêne - BP 50072 - HERBLAY - 95612 CERGY PONTOISE CEDEX - FRANCE