



STARLIGHT PRO

■ stop too complex

Many polymerization units are difficult to adjust and complicated in their handling – which makes it really hard to integrate them into your workflow.

► start to work

Well, what about just starting to work? Due to its cordless design, you can use the lightweight starlight pro wherever you want to – and its one-button-handling couldn't be any easier.



► lithium-ion battery

The lithium-ion battery provides you with 320 working cycles of 10 seconds - after only 90 minutes of charging.



► click-clack connection

Exclusive connection in metal ensures an optimal light beam transmission.



► lightweight

Only 105 g.



► high performance level

The power of more than 1.400 mW/cm² lets you harden a 2-mm thick layer in only 10 seconds.

► optical fiber

The sterilizable optical fiber (8 mm) is made of many single threads for maximum light power.



► light spectrum

The light spectrum from 440 to 480 nm is focused on camphorquinone, a photoinitiator used in 98 % of dental materials.

